

# Freedom of choice

We offer a full line of oxygen concentrators

## EverFlo Q ultra quiet oxygen concentrator with all of the features of our EverFlo, plus:

- Provides quiet sound level of 40 dBA (under typical use)
- Reduces noise complaints and associated costs and time
- Differentiates your oxygen business to physicians



## SimplyGo and SimplyGo Mini Portable Oxygen Concentrators

- SimplyGo offers 6 pulse settings and 2LPM continuous flow oxygen in 10 pounds
- SimplyGo Mini offers 5 pulse settings in 5 pounds
- Both SimplyGo and SimplyGo Mini have easy to use non-medical designs that plug into any AC or DC power source
- SimplyGo and SimplyGo Mini comply with all FAA acceptance criteria for POCs



Complies with RTCA/DO-160G section 21, category M  
The manufacturer of this product has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft.



PHILIPS  
RESPIRONICS

EverFlo

Home oxygen system



# Oxygen simplified

© 2017 Koninklijke Philips N.V. All rights reserved.  
Specifications are subject to change without notice.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Caution: US federal law restricts these devices to sale by or on the order of a physician.

1010 Murry Ridge Lane, Murrysville, PA 15668  
800-345-6443 • 724-387-4000

[www.philips.com/respironics](http://www.philips.com/respironics)

KR 05/22/17 MCI 4107943 PN 1038667



# Ideal. Oxygen. Concentrators.

For you and your patients.

## EverFlo concentrators can reduce the cost of oxygen therapy.

Philips Respironics recognizes the growing need for cost-effective and reliable oxygen concentrators that can help you manage the cost of providing oxygen therapy. That's why we have EverFlo, a series of forward thinking oxygen products designed to simplify your oxygen business.

### Expect industry-leading quality and reliability

EverFlo stationary oxygen concentrators can save you time and money over the course of ownership because of their design and quality. Field proven since 2007, there have been more than 2.5 million EverFlo units sold worldwide.\*

Assembled in the United States and backed by our deep design knowledge, EverFlo is developed using only high-quality parts. With a recessed flow meter, metal cannula connector, strong plastic casing, proven sieve bed, and Pressure Swing Adsorption (PSA) technology, you may notice a reduction in repairs, complaints, and unscheduled visits.

### Enjoy ease of service

EverFlo was designed keeping maintenance costs in mind. For example, EverFlo has just two filters. The external filter can be changed in minutes without special skills or tools and the internal micro disk filter is simply checked during routine maintenance procedures.

### Differentiate your business

EverFlo concentrators are backed by value-added tools. From in-person and online training to financing options, you have what you need to succeed. And, since Philips Respironics offers a wide range of sleep and home respiratory products and comprehensive support, we can be the source for many of your product needs.

\*Data on file Philips Respironics 2017

### Easy to store, transport, and deliver.

EverFlo is compact, has a streamlined and integrated handle, and weighs only 31 pounds. Less weight means EverFlo is easier to move around and can translate into lower shipping costs.

### Give your patients what they want.

At just 31 pounds, your patients will appreciate the compact, non-medical design of the EverFlo oxygen concentrator because it will fit right into their home environment. The unit is designed to be small, light, easy to transport and will fit into the trunk of a car. EverFlo has an easy-to-mount humidifier bottle connection and no patient filter.



**Metal cannula**  
Connector for greater durability

**Molded handle**  
that's easy to pick up and transport with one or both hands

**Interface panel**  
with on/off switch and alarm indicators

**Recessed flow meter**  
reduces the risk of breakage



**Internal tubing storage**  
keeps tubing conveniently available for use



**Humidifier bottle holder**  
compatible with all bottle styles

### EverFlo product specifications

Input voltage	120 VAC +/- 10%
Input frequency	60 hz
Average power consumption	350 watts
Oxygen concentration* (at 5 lpm)	93% (+/- 3%)
Liter flow	0.5 to 5 liters per minute
Weight	31 lbs (14 kg)
Sound level	45 dBA typical
Dimensions	23" (584 mm) h, 15" (381 mm) w, 9.5" (241 mm) d
Required maintenance	Compressor intake filter: change every 2 years
OPI (oxygen percentage indicator) alarm levels	Low oxygen: 82%, very low oxygen: 70%
Operating temperature	55°F to 90°F (12°C to 32°C)
Storage/transport humidity	-30°F to 160°F (-34°C to 71°C) up to 95% relative humidity
Operating humidity	Up to 95% relative humidity
Outlet pressure	5.5 PSI
Operating altitude	0 to 7,500' (0 to 2,286 m)

### EverFlo ordering information

Concentrator systems	Part number
EverFlo oxygen concentrator	1020000
EverFlo oxygen concentrator with OPI (oxygen percentage indicator)	1020001
EverFlo Q oxygen concentrator	1020015
EverFlo Q oxygen concentrator with OPI (oxygen percentage indicator)	1020014

Concentrator accessories	Part number
Compressor intake filter	1038831
Micro disk filter	520
Humidifier connector tube	1039642
Stand-alone low-flow flow meter	H537-1
Integrated low-flow flow meter	H644
Flow meter locking kit	528
EverFlo manual in English, French, Spanish, and Portuguese	1020121
EverFlo service manual	1039055

Extended warranty	Part number
Two-year extended warranty with compressor	1038571

### Pair EverFlo with SimplyGo or SimplyGo Mini portable oxygen concentrators for a complete oxygen solution.

When SimplyGo or SimplyGo Mini is teamed up with EverFlo, you can create an oxygen-generating system that will meet the needs of nearly all on-the-go patients while helping to reduce your cost of providing oxygen therapy.

